

SpermMorph

"Strict Morphology" Pre-Stained Human Sperm
Quality Control Microscope Slides for Manual Application

QC Protocol

Equipment Needed: 1000x light microscope with oil and micrometer.

Set Includes: Slide 1: Fertile; > 14% normal forms (normal pattern)
Slide 2: Subfertile; 5 - 14% normal forms (g-pattern)
Slide 3: Subfertile; 0 - 4% normal forms, proven failed fertilization (pattern)

Per request, non-labeled slides may be purchased for blind surveys.

Protocol: Performing QC in your laboratory is extremely important because it improves precision and accuracy which translates into improved clinical results." A minimal 100 cells (100-200) are required for complete analysis with results expressed using "Tygerberg Strict Criteria." Low variability per slide exists with inter-technician (0.86 c.c) and intra-technician (0.96 c.c.) using the "Spearman Rank Correlation Coefficient."

(c.c. = Correlation Coefficient)

Strict Criteria for a Normal Sperm Cell^{2,5}

Length: 4 - 5.5 microns⁶
Width: 2.5 - 3.5 microns⁶
Shape: Smooth Oval (acorn)^{2,5}
Acrosome: 40 - 70% of head^{2,5}
Midpiece: 6 - 10 microns
Tail: +/- 45 microns
Other Criteria: Slightly amorphous (borderline) forms = abnormal^{2,5,6}
No severe neck, midpiece or tail defects are allowed^{2,5,6}

¹ Oehninger S, Acosta AA, Morshedi M, veeck LL, Swanson RJ, Simmons K, Rosenwaks Z: Corrective Measures and pregnancy outcome in in vitro fertilization in patients with severe sperm morphology abnormalities. *Fertil Steril* 1988;50:283

² Kruger TF, Menkveld R, Stander FSH, Lombard CJ, Van der Merwe JP, Van Zyl JA, Amith K: Sperm Morphology features as prognostic factor in in vitro fertilization. *Fertil Steril* 1986; 46:1118.

³ Van der Merwe JP, Kruger TF, Swart Y, Lombard CJ: The role of oocyte maturity in the treatment of infertility because of teratozoospermia and normozoospermia with gamete intrafallopian transfer. *Fertil Steril* 1992; Vol 58, No 3.

⁴ Kruger TF, Swanson RJ, Hamilton M, Simmons K, Acosta AA, Matta JF, Oehninger S, Morshedi M: Abnormal sperm morphology and other semen parameters related to the results of the hamster oocyte human sperm penetration assay. *Int J Androl* 1988;11:107.

⁵ Menkveld R, Stander FSH, Kotze TJvW, Kruger TF, Van Zyl JA: The evaluation of morphological characteristics of human spermatozoa according to stricter criteria.

⁶ World Health Organization. Laboratory manual for the examination of human semen and sperm-cervical mucus interaction. Cambridge University Press, 3rd Edition 1992.

Manufactured by MQ Medical, Tygerberg-Capetown, South Africa