



## Genitor™

### Description

Double-lined transparent polyethylene cannula:

- 17 cm long
- OD: 1.66 mm
- ID: 0.65 mm
- protruding inner lining forms a very flexible blunt tip with two side ports
- standardized adaptor allows connection to a syringe

### Indications: Intrauterine Artificial Insemination

The narrow diameter of the Genitor™ reduces the residual (dead) volume by 80% compared to the standard model. This is particularly useful for donor sperm or when a very small quantity is available for insemination. The position of the apertures allows deposition of the semen close to the tubal orifices.

### CCD Code

120400M Genitor™ (box of 25 units)

[www.sepalreprodevices.com](http://www.sepalreprodevices.com)

www

410 Boston Post Road, Suite 21  
Sudbury, Massachusetts 01776 USA  
Email: [info@sepalreprodevices.com](mailto:info@sepalreprodevices.com)  
Main: 978.579.9389  
Fax: 978.579.9407  
Toll Free: 877.546.4223



Manufactured by Laboratoire C.C.D., Paris, France