



Intrauterine Insemination Cannula with mandrel

Description

Double-lined transparent polyethylene cannula:

- 17 cm long
- OD: 2.1 mm
- ID: 1.7 mm
- protruding inner lining forms a very flexible blunt tip with two side ports
- standardized adaptor allows connection to a syringe
- the cannula is provided with a malleable steel mandrel covered with a polyethylene sheath

Indications: Intrauterine Artificial Insemination

This cannula is adapted for A.I. through poorly accessible and tight cervixes. Insertion of the cannula is easier because of the mandrel. The cannula allows for easy and atraumatic access to the fundus of the uterus. The position of the apertures allows deposition of the semen close to the tubal orifices.

CCD Code

12040AM Intrauterine Insemination Cannula with Mandrel (box of 25 units)

www.sepalreprodevices.com

410 Boston Post Road, Suite 21
Sudbury, Massachusetts 01776 USA
Email: info@sepalreprodevices.com
Main: 978.579.9389
Fax: 978.579.9407
Toll Free: 877.546.4223



Manufactured by Laboratoire C.C.D., Paris, France