



# Standard Intrauterine Insemination Cannula

## Description

Double-lined transparent polyethylene cannula:

- 17 cm long
- OD: 2.1 mm
- ID: 1.7 mm
- protruding inner lining forms a very flexible blunt tip with two side ports
- standardized adaptor allows connection to a syringe

## Indications: Intrauterine Artificial Insemination

Although the outer diameter of this catheter offers a small degree of rigidity in cannulating the cervix, its flexibility facilitates atraumatic access to the uterine fundus. The position of the apertures allows deposition of the semen close to the tubal orifices.

## CCD Code

120400S Standard Intrauterine Insemination Cannula (box of 25 units)

[www.sepalreprodevices.com](http://www.sepalreprodevices.com)

www

410 Boston Post Road, Suite 21  
Sudbury, Massachusetts 01776 USA  
Email: [info@sepalreprodevices.com](mailto:info@sepalreprodevices.com)  
Main: 978.579.9389  
Fax: 978.579.9407  
Toll Free: 877.546.4223

CE 0120 STERILE R Mouse Embryo toxicity tested MEA Endo Toxine Level tested LAL

Manufactured by Laboratoire C.C.D., Paris, France