



FIV 2

FIV 4


FIV 6

Laboratory Countertop Incubators

*Designed to Optimize In Vitro Fertilization and Reproductive Biology
Improved accuracy, adjustability, and consistency during incubation*

Features

- compact size frees counter space for other uses
- constant temperature control facilitated by display counter
- adjustable 37°C temperature checkpoint
- accuracy to $\pm 0.2^{\circ}\text{C}$ provided by temperature detector
- constant control for gas content
- heated covers
- thermostatic humidifier



Order Codes

CA-FIV-2	Laboratory Countertop Incubator
CA-FIV-4	Laboratory Countertop Incubator
CA-FIV-6	Laboratory Countertop Incubator