



Frydman™ FEP Ultrasoft Echo Metal

Description

The set includes 2 catheters:


- Catheter N° 1: Fluoro Ethylene Propylene Introducer:
 - 17 cm long
 - ID: 1.90 mm
 - OD: 2.30 mm, with positioning markings at 1, 2, 3, 4, and 5 cm from its distal end
- Catheter N° 2: Polyurethane Re-Implantation Catheter:
 - 25 cm long with a soft rounded tip
 - ID: 0.70 mm
 - OD: 1.55 mm
 - the distance between the end of the introducer and the tip of the catheter is 5 cm
 - with 5 markings spaced 1 cm apart in its lower part
 - the small pharmaceutical grade stainless steel ring has sandblasted corners to ensure a smooth surface

Indications: IVF: Embryo Transfer

The Fluoro Ethylene Propylene (FEP) material provides extraordinary memory that can be shaped prior to insertion and will maintain its integrity. The partial metallic reinforced inner catheter is designed for easier transfers.

C.C.D. Code

1321501 Frydman™ FEP Ultrasoft Echo Metal (box of 10 units)



www.sepalreprodevices.com

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA

main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223

email: info@sepalreprodevices.com