



# Frydman™ FEP Ultrasoft

*with or without Echogenic Tip*

## Description

*The set includes 2 catheters:*

- Catheter N° 1: Fluoro Ethylene Propylene (FEP) Introducer:
  - 17 cm long
  - ID: 1,90 mm
  - OD: 2,30 mm, with positioning markings at 1, 2, 3, 4, and 5 cm from its distal end
- Catheter N° 2: Polyurethane Re-Implantation Catheter:
  - 25 cm long with a soft rounded tip
  - ID: 0,70 mm
  - OD: 1,55 mm
- the distance between the end of the introducer and the tip of the catheter is 5 cm
- with 5 markings spaced 1 cm apart in its lower part
- the small pharmaceutical grade stainless steel ring has sandblasted corners to ensure a smooth surface
- Available with a separate FEP Introducer with Guide
- Available in 20 cm length


## Indications: IVF: Embryo Transfer

The Fluoro Ethylene Propylene material provides extraordinary memory that can be shaped prior to insertion and will maintain its integrity, while the soft inner catheter minimizes trauma to the cervical canal and endometrium.

## C.C.D. Codes

1324501	Frydman™ FEP Ultrasoft 25 cm (box of 10 units)
1322301	Frydman™ FEP Ultrasoft Echo 25 cm (box of 10 units)
1321301	Frydman™ FEP Ultrasoft 20 cm (box of 10 units)
1322701	Frydman™ FEP Introducer with Guide 25 cm (box of 10 units)
1350301	Frydman™ FEP Introducer with Guide 20 cm (box of 10 units)

*Manufactured by Laboratoire C.C.D., Paris, France*



[www.sepalreprodevices.com](http://www.sepalreprodevices.com)

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA  
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223  
email: [info@sepalreprodevices.com](mailto:info@sepalreprodevices.com)