



Frydman™ FEP Ultrasoft Mock

with or without Echogenic Tip

Description

The set includes 2 catheters:

- Catheter N° 1: Fluoro Ethylene Propylene (FEP) Introducer:
 - 17 cm long
 - ID: 1,90 mm
 - OD: 2,30 mm, with positioning markings at 1, 2, 3, 4, and 5 cm from its distal end
- Catheter N° 2: Polyurethane Re-Implantation Catheter:
 - 25 cm long with a soft rounded tip
 - ID: 0,70 mm
 - OD: 1,55 mm
- the distance between the end of the introducer and the tip of the catheter is 5 cm
- with 5 markings spaced 1 cm apart in its lower part
- the small pharmaceutical grade stainless steel ring has sandblasted corners to ensure a smooth surface
- Available with a separate FEP Introducer with Guide
- Available in 20 cm length

Indications: IVF: Embryo Transfer

With the same format and softness as the Ultrasoft FEP Frydman™, but sealed to their extremities, the Ultrasoft Mock FEP enables evaluation of the permeability of the cervix, its front or backward orientation, and the possible presence of cervical obstructions. The FEP Ultrasoft Mock accurately reproduces real transfer conditions without cervical dilation since its diameter is identical to the Frydman™ FEP Ultrasoft.

C.C.D. Codes

1324601	Frydman™ FEP Ultrasoft Mock 25 cm (box of 10 units)
1322401	Frydman™ FEP Ultrasoft Echo Mock 25 cm (box of 10 units)
1323001	Frydman™ FEP Ultrasoft Mock 20 cm (box of 10 units)
1322701	Frydman™ FEP Introducer with Guide 25 cm (box of 10 units)
1350301	Frydman™ FEP Introducer with Guide 20 cm (box of 10 units)

Manufactured by Laboratoire C.C.D., Paris, France



www.sepalreprodevices.com

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223
email: info@sepalreprodevices.com