



T.D.T.™ Set (Tight Difficult Transfer)

Description

The set includes 2 catheters:

- Catheter N° 1: Polypropylene Introducer:
 - 17,5 cm long
 - flexible tip, 4,5 cm long
 - two hysterometric markings positioned at 5,5 and 6,5 cm from its distal end
 - provided with a plastic coated metal mandrel
- Catheter N° 2: Ultra-thin Polyethylene Re-Implantation Catheter:
 - 23,5 cm long
 - on a steel micro-tube
 - OD: 1,0 mm
 - ID: 0,6 mm
 - with a positioning mark at 2 cm from its lower part


Indications: IVF: Embryo Transfer

The T.D.T.™ (Tight Difficult Transfer) set reduces the risk of mucus accumulation in the transfer catheter that would impede the expulsion of the embryo.

CCD Code

1308000 T.D.T.™ Set (box of 10 units)

Manufactured by Laboratoire C.C.D., Paris, France



www.sepalreprodevices.com

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223
email: info@sepalreprodevices.com