



Luer Lock Needle for Ovum Recovery

Description

- Single flow needle, 17 gauge O.D., 18 gauge I.D. with a double-cutting beveled point and an echogenic end part on 2 cm, and over 5 mm from the distal 3rd cm.
- Available in 2 lengths: 25 cm, 30 cm.
- Also available in a 16 gauge, 35 cm length, ID: 1.30 mm, OD: 1.60 mm.


Indications: IVF: Embryo Transfer

Collection of oocytes using celioscopic approach or ultrasound guidance (transvesical, transvaginal, periurethral or endovaginal).

C.C.D. Codes

1403002	Luer Lock Needle, 17 gauge, 25 cm (box of 10 units)
1403001	Luer Lock Needle, 17 gauge, 30 cm (box of 10 units)
1311101	Luer Lock Needle, 16 gauge, 35 cm (box of 10 units)

Manufactured by Laboratoire C.C.D., Paris, France



www.sepalreprodevices.com

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223
email: info@sepalreprodevices.com