



Short Frydman™ Set and Long Frydman™ Set

Short Frydman™ Set Description

The set includes 2 catheters:

- Catheter N° 1: Polypropylene Introducer:
 - 14 cm long
 - OD: 2.2 mm with positioning markings at 1, 2, 3, 4, 5, 6, 7 cm from its distal end
 - with a green mandrel and a movable sliding stopper
- Catheter N° 2: Polyethylene Re-Implantation Catheter:
 - 21 cm long
 - ID: 1.1 mm
 - OD: 1.6 mm
 - with a steel tube and three markings spaced 1 cm apart in its lower part

Long Frydman™ Set Description

The set includes 2 catheters:

- Catheter N° 1: Polypropylene Introducer:
 - 20 cm long
 - OD: 2.2 mm with positioning markings at 1, 2, 3, 4, 5, 6, 7 cm from its distal end
 - with a green mandrel and a movable sliding stopper
- Catheter N° 2: Polyethylene Re-Implantation Catheter:
 - 25.5 cm long
 - ID: 1.1 mm
 - OD: 1.6 mm
 - with a steel tube and one marking on its lower part


Indications: IVF: Difficult Embryo Transfer

This set may be used in the presence of a narrow or curved cervix that makes it difficult to introduce a standard Frydman™ catheter.

C.C.D. Codes

1307000 Short Frydman™ Set (box of 10 units)
1303000 Long Frydman™ Set (box of 10 units)

Manufactured by Laboratoire C.C.D., Paris, France



www.sepalreprodevices.com

410 Boston Post Road • Suite 21 • Sudbury, Massachusetts • 01776 • USA
main: 978.579.9389 • fax: 978.579.9407 • toll free: 877.546.4223
email: info@sepalreprodevices.com