



Soft Frydman™ Set

*with or without Mandrel and/or Echogenic Tip
also available as a Mock Catheter with or without Echogenic Tip
also available, separate Introducer with Guide*

Description

The set includes 2 catheters:

- Catheter N° 1: Polypropylene Introducer:
 - 14,5 cm long
 - OD: 2,2 mm with positioning markings at 1, 2, 3, 4, 5, 6, 7 cm from its distal end
- Catheter N° 2: Polyethylene Re-Implantation Catheter:
 - 23 cm long
 - ID: 1,53 mm
 - OD: 0,7 mm
 - with five markings spaced 1 cm apart in its lower part


Indications: IVF: Embryo Transfer

The model with mandrel may be used in the presence of a narrow or curved cervix that would make it difficult to introduce the model without mandrel.

C.C.D. Codes

| | |
|---------|--|
| 1307100 | Frydman™ Soft with Mandrel (box of 10 units) |
| 13070S0 | Frydman™ Soft without Mandrel (box of 10 units) |
| 1323401 | Frydman™ Soft Echo with Mandrel (box of 10 units) |
| 1323301 | Frydman™ Soft Echo without Mandrel (box of 10 units) |
| 1323601 | Frydman™ Soft Mock (box of 10 units) |
| 1323901 | Frydman™ Soft Echo Mock (box of 10 units) |
| 1321101 | Frydman™ Polypropylene Introducer with Guide (box of 20 units) |

Manufactured by Laboratoire C.C.D., Paris, France



www.sepalreprodevices.com

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223
email: info@sepalreprodevices.com