



# Ultrasoft Frydman™ Set

*with or without Mandrel and/or Echogenic Tip  
also available as a Mock Catheter with or without Echogenic Tip  
also available, separate Introducer with Guide*

## Description

*The set includes 2 catheters:*

- Catheter N° 1: Polypropylene Introducer:
  - 14,5 cm long
  - OD: 2.2 mm with positioning markings at 1, 2, 3, 4, 5, 6, 7 cm from its distal end
- Catheter N° 2: Polyurethane Re-Implantation Catheter:
  - 23 cm long
  - ID: 1,53 mm
  - OD: 0,7 mm
  - With five markings spaced 1 cm apart in its lower part


## Indications: IVF: Embryo Transfer

The model with mandrel may be used in the presence of a narrow or curved cervix that would make it difficult to introduce the model without mandrel.

## C.C.D. Codes

131230301	Frydman™ Ultrasoft with Mandrel (box of 10 units)
131220301	Frydman™ Ultrasoft without Mandrel (box of 10 units)
1322801	Frydman™ Ultrasoft Echo with Mandrel (box of 10 units)
1323501	Frydman™ Ultrasoft Echo without Mandrel (box of 10 units)
1323701	Frydman™ Ultrasoft Mock (box of 10 units)
1323801	Frydman™ Ultrasoft Echo Mock (box of 10 units)
1321101	Frydman™ Polypropylene Introducer with Guide (box of 20 units)

*Manufactured by Laboratoire C.C.D., Paris, France*



[www.sepalreprodevices.com](http://www.sepalreprodevices.com)

410 Boston Post Road · Suite 21 · Sudbury, Massachusetts · 01776 · USA  
main: 978.579.9389 · fax: 978.579.9407 · toll free: 877.546.4223  
email: [info@sepalreprodevices.com](mailto:info@sepalreprodevices.com)